


Olympus Upright Microscope Brightfield


1, Turn the equipment on

- a. Power switch for microscope is on the right at the back
- b. Adjust the light with the slider (move to the camera icon)




2. Perform Kohler Illumination (*for step-by-step pictorial view see following page*)

- a. Put a slide on the stage. Use the 10x objective.
- b. Open the condenser and the field diaphragm all the way.
- c. Move the condenser up close to the slide, using the condenser focus.
- d. Focus the sample. Use the neutral density filters if it's too bright.
- e. Close the field diaphragm as far as it goes. A small circle of light should be seen in the eyepieces.
- f. Center the small circle using the condenser centering rods.
- g. Focus it until the blade edges can be seen using the condenser focus. It will probably have to be re-centered again after this.
- h. Open the field diaphragm all the way.
- i. Take off one of the eyepieces and close the condenser by about 10%.



3. Running the Scion VisCapture software


- a. Open the Scion VisCapture software
- b. **Image > Show Properties** to pull up the **Camera Properties** (this can be left up)
- c. Select the **Brightfield** profile under **Configuration Select**
- d. Perform White Balance
 - i. Click on the **Color Balance** tab
 - ii. Select an area off the sample
 - iii. Click **AWB**
- e. Under **Settings**, adjust the Gain and Exposure to lighten or darken
- f. Set **Gamma** to about 0.6 (under **Color** tab)
- g. Click on the **Multi-Frame** tab and set **Average** to **4**
- h. Taking a picture
 - i. **Click** on the image or go **Image > Snap**
 - ii. **File > Save**
 - iii. **Image > Start Live Capturing** or press **Ctrl +G** to go back to live mode




Olympus Upright Microscope Fluorescence

1, Turn the equipment on


- a. Power switch for microscope is on the right back
- b. Turn the lamp on (box behind microscope)
- c. Open the shutter by pushing it **IN**
- d. Select cube






2. Running the Scion VisCapture software

- a. Open the Scion VisCapture software
- b. **Image > Show Properties** to pull up the **Camera Properties** (this can be left up)
- c. Select a **Fluorescence** profile under **Configuration Select**
- d. Under **Settings**, adjust the Gain and Exposure to lighten or darken
- e. Set **Gamma** to about 1 (under **Color** tab)
- f. Check the color gain
 - i. **Green Fluorescence**
Green Gain at 10.2 db, Red and Blue at -2.0 db
 - ii. **Red Fluorescence**
Red Gain at 10.2 db, Green and Blue at -2.0 db
 - iii. **Blue Fluorescence**
Blue Gain at 10.2 db, Red and Green at -2.0 db
- f. Click on the **Multi-Frame** tab and set **Average** to 4
- g. Taking a picture
 - i. **Click** on the image or go **Image > Snap**
 - ii. **File > Save**
 - iii. **Image > Start Live Caputring** or press **Ctrl +G** to go back to live mode



Default Settings for Scion VisCapture Software



Brightfield Settings

Gamma at 0.6

Red, Green, and Blue Gain in a rough line

Fluorescence Settings

Gamma at 1

Fast Fluorescence

Bin makes smaller images but is faster (used for dim pictures)

Green Fluorescence

Green Gain at 10.2 db, Red and Blue at -2.0 db

Red Fluorescence

Red Gain at 10.2 db, Green and Blue at -2.0 db

Only Red Boost checked

Blue Fluorescence

Blue Gain at 10.2 db, Red and Green at -2.0 db

Only Blue Boost checked